

---

# LabCorp OTS Web Services completeCollection2

---

Revision: 1.0

Date: 1/25/2016



**Occupational Testing Services**

**Information Systems  
completeCollection2**

**Revision History**

<b>Revision</b>	<b>Date</b>	<b>Author</b>	<b>Description</b>
1.0	01/25/2016	Piotr Gucwa	Initial Draft

**System:** OTS Web Services

**Title:** completeCollection2

**Filename:** OTS\_WS\_ completeCollection2.doc

**Revision:** 1.0

**Date:** 1/25/16

**Page:** 2 of 3

```
class CollectionRequest {
    xsd:string(10) specimenID;
    xsd:string(2) specimenType;
    //U-urine, H-hair, O-oral fluid, BU-blood and urine
    xsd:string(10) collectionSiteID;
    xsd:dateTime collectionCompletedDateTime;
    xsd:string(25) collectorNameFirst;
    xsd:string(1) collectorNameMiddleInitial;
    xsd:string(35) collectorNameLast;
    xsd:string(1) donorIdentifiedby;
    //G-photo id, E-employer representative
    xsd:boolean donorRefusedtoTest;
    xsd:boolean nospecimenprovided;
    xsd:boolean withinTemperatureRange;
    //applicable to urine specimen
    xsd:int temperature;
    //deg Farenheit - specimen temperature
    xsd:boolean collectedFromHead;
    xsd:boolean splitSpecimenCollected;
    xsd:boolean collectionObserved;
    xsd:string(2) POCTResults;
    //NG-negative, NN-non-negative, NR-no result
    xsd:string(8) deviceLotNumber;
    xsd:dateTime deviceExpirationDate;
    xsd:string(2) deliveryService;
    //LC-Labcorp Courier, FX-fedex, OT-other
    xsd:string(200) collectorRemarks;
}
```

/\*  
This call is used by non Labcorp collection sites to pass collection information to Labcorp after collection was completed.  
When the specimen arrives at lab we will scan the specimen number and populate the specimen data with the data passed in this call.  
Before calling completeCollection2(), donor should be registered in Labcorp system. Labcorp registration number has to be passed to this call.  
\*/

```
void completeCollection2(
    xsd:string userId,
    xsd:string password,
    int labcorpRegistrationNumber,
    CollectionRequest collection
) throws WsException;
```